

CLIMBUP® Insect Interceptor is available in 4 sizes:

DIMENSIONS (inches)

| | Shape | Height | Inner well | Total Size |
|----------------|--------------|---------------|-------------------|-------------------|
| CLIMBUP | round | 1" | 4.25" | 7" |

This is the original model that will hold most bed legs and is low enough to go under furniture skirts.

| | | | | |
|------------------|-----------|----|--------------|---------------|
| CLIMBUPXL | rectangle | 2" | 4.62"X 6.62" | 7.36" X 9.35" |
|------------------|-----------|----|--------------|---------------|

This is an extra large model to hold the legs of attached headboards and bed legs within the same receptacle.

| | | | | |
|--------------------|-------|----|------|------|
| CLIMBUPtank | round | 2" | 5.5" | 9.5" |
|--------------------|-------|----|------|------|

This is a model which will hold a 5 lb. compressed carbon dioxide tank for slow release of carbon dioxide when human breathe is not present to act as the attractant into the CLIMBUP.

| | | | | |
|-------------------|-------|----|----|-------|
| CLIMBUPcup | round | 2" | 1" | 4.75" |
|-------------------|-------|----|----|-------|

The inner well of this model will hold very small cot legs. It may also be used adjacent to bed where someone is sleeping to monitor activated bed bugs but it will not serve as a barrier in this application.